
Product Specifications

LBL-Dye A715 NHS Ester

<https://proimaging.fr/home/39-lbl-dye-a715-nhs-ester-0210704.html>

Ordering Information

Product Number: 0210704

Product Name: LBL-Dye A715 NHS Ester

Product Reference: BM105

Unit Size: 1mg / 5 mg

Storage Conditions: Freeze (≤ -15 °C), Desiccated, Avoid Light

Expiration Date: 12 months upon receiving

Chemical and Spectral Properties

Molecular Weight: ~ 1212

Solubility:


-PBS Buffer (pH 7.4) + 1% polysorbate 80

-DMSO

-Water

Max Absorption: 690 nm

Max Emission: 715 nm



Analysis medium : PBS buffer (pH 7.4) + 1% polysorbate 80

Application Notes

LBL-Dye A715 NHS Ester is an activated bright dye with high photostability.

This dye is water-soluble and ideal for labeling sensitive proteins, peptides, and antibodies.

Product for Research Use Only. Not for Human Use.