

---

## Product Specifications

### LBL-Dye A715 NHS Ester

<https://proimaging.fr/product/lbl-dye-a715-nhs-ester/>

### Product Information

**Product Name:** LBL-Dye A715 NHS Ester

**Product Reference:** BM105

**Unit Size:** 1mg

**Shipping Conditions:** Room temperature for up to 3 weeks

**Storage Conditions:** Store 18 months at -20°C and protect from light

**Molecular Weight:** ~ 1212

**Max Absorption:** 690 nm


**Max Emission:** 715 nm

**Coupling medium:**

pH 8.5-9.0 phosphate buffer  
+ 10% DMF (or DMSO)

**Solubility:**

-PBS Buffer (pH 7.4) + 1% polysorbate 80  
-DMSO  
-Water



### Product Description

LBL-Dye A715 NHS Ester is a highly photostable, activated bright dye.

This water-soluble dye is ideal for labeling sensitive proteins, peptides, and antibodies.

**Product for Research Use Only. Not for Human Use.**