

Product Specifications

LBL-Dye A820 NHS Ester

<https://proimaging.fr/home/40-lbl-dye-a820-nhs-ester-0210705.html>

Product Information

Product Name: LBL-Dye A820 NHS Ester

Product Reference: BM146

Unit Size: 1mg

Shipping Conditions: Room temperature for up to 3 weeks

Storage Conditions: Store 18 months at -20°C and protect from light

Molecular Weight: ~ 1244

Max Absorption: 800 nm


Max Emission: 820 nm

Coupling medium:

pH 8.5-9.0 phosphate buffer
+ 10% DMF (or DMSO)

Solubility:

-PBS Buffer (pH 7.4) + 1% polysorbate 80
-DMSO
-Water



Product Description

LBL-Dye A820 NHS Ester is a highly photostable, activated bright dye.

This water-soluble dye is ideal for labeling sensitive proteins, peptides, and antibodies.

Product for Research Use Only. Not for Human Use.